



EXTERNAL HEXAGON

4.1 REGULAR CONNECTION

CYLINDRICAL SHAPE



EC




IMPIANTO
CON MOUNTER
WITH MOUNTER


EC 34

CYLINDRIC
external hexagon - full treatment


NARROW CYLINDRICAL SHAPE



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3408	Ø 3.4	8.5 mm	Ø 4.1	Ø 2.5
EC 3410	Ø 3.4	10 mm	Ø 4.1	Ø 2.5
EC 3411	Ø 3.4	11.5 mm	Ø 4.1	Ø 2.5
EC 3413	Ø 3.4	13 mm	Ø 4.1	Ø 2.5




REGULAR CYLINDRICAL SHAPE



- 0.5 mm Machined
- Cylindric
- Screw Pitch 1.2 mm DT
- Tapered 3 mm




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3808	Ø 3.8	8.5 mm	Ø 4.1	Ø 2.7
EC 3810	Ø 3.8	10 mm	Ø 4.1	Ø 2.7
EC 3811	Ø 3.8	11.5 mm	Ø 4.1	Ø 2.7
EC 3813	Ø 3.8	13 mm	Ø 4.1	Ø 2.7
EC 3815	Ø 3.8	15 mm	Ø 4.1	Ø 2.7




EC 42

CYLINDRIC
external hexagon - full treatment

REGULAR CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 4206	Ø 4.2	6,0 mm	Ø 4.1	Ø 3.1
EC 4208	Ø 4.2	8.5 mm	Ø 4.1	Ø 3.1
EC 4210	Ø 4.2	10 mm	Ø 4.1	Ø 3.1
EC 4211	Ø 4.2	11.5 mm	Ø 4.1	Ø 3.1
EC 4213	Ø 4.2	13 mm	Ø 4.1	Ø 3.1



CYLINDRIC
external hexagon - full treatment

EC 50

WIDE CYLINDRICAL SHAPE



0.7 mm
Machined

Cylindric

Screw Pitch
1.6 mm DT


Tapered
3 mm



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 5006	Ø 5.0	6.0 mm	Switch Ø 4.1	Ø 3.9
EC 5008	Ø 5.0	8.5 mm	Switch Ø 4.1	Ø 3.9
EC 5010	Ø 5.0	10 mm	Switch Ø 4.1	Ø 3.9
EC 5011	Ø 5.0	11.5 mm	Switch Ø 4.1	Ø 3.9
EC 5013	Ø 5.0	13 mm	Switch Ø 4.1	Ø 3.9

SOFT AND STD BONE - D2/D3/D4


HARD BONE - D1




EC 34 HT

CYLINDRIC HT
external hexagon - half treatment

NARROW CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3408 HT	Ø 3.4	8.5 mm	Ø 4.1	Ø 2.5
EC 3410 HT	Ø 3.4	10 mm	Ø 4.1	Ø 2.5
EC 3411 HT	Ø 3.4	11.5 mm	Ø 4.1	Ø 2.5
EC 3413 HT	Ø 3.4	13 mm	Ø 4.1	Ø 2.5




CYLINDRIC HT
external hexagon - half treatment

EC 38 HT

REGULAR CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3808 HT	Ø 3.8	8.5 mm	Ø 4.1	Ø 2.7
EC 3810 HT	Ø 3.8	10 mm	Ø 4.1	Ø 2.7
EC 3811 HT	Ø 3.8	11.5 mm	Ø 4.1	Ø 2.7
EC 3813 HT	Ø 3.8	13 mm	Ø 4.1	Ø 2.7
EC 3815 HT	Ø 3.8	15 mm	Ø 4.1	Ø 2.7




EC 42 HT

CYLINDRIC HT
external hexagon - half treatment

REGULAR CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 4208 HT	Ø 4.2	8.5 mm	Ø 4.1	Ø 3.1
EC 4210 HT	Ø 4.2	10 mm	Ø 4.1	Ø 3.1
EC 4211 HT	Ø 4.2	11.5 mm	Ø 4.1	Ø 3.1
EC 4213 HT	Ø 4.2	13 mm	Ø 4.1	Ø 3.1



CYLINDRIC HT
external hexagon - half treatment

EC 50 HT

WIDE CYLINDRICAL SHAPE



3.5 mm
Machined

Cylindric

Screw Pitch
1.4 mm DT

Tapered
3 mm



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 5008 HT	Ø 5.0	8.5 mm	Switch Ø 4.1	Ø 3.9
EC 5010 HT	Ø 5.0	10 mm	Switch Ø 4.1	Ø 3.9
EC 5011 HT	Ø 5.0	11.5 mm	Switch Ø 4.1	Ø 3.9
EC 5013 HT	Ø 5.0	13 mm	Switch Ø 4.1	Ø 3.9

SOFT AND STD BONE - D2/D3/D4

HARD BONE - D1

