



EXTERNAL HEXAGON

4.1 REGULAR CONNECTION

TAPERED
SHAPE




EK




IMPIANTO
CON MOUNTER
WITH MOUNTER



REGULAR TAPERED SHAPE



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EK 3808	Ø 3.8	8.5 mm	Ø 4.1	Ø 1.8
EK 3810	Ø 3.8	10 mm	Ø 4.1	Ø 1.8
EK 3811	Ø 3.8	11.5 mm	Ø 4.1	Ø 1.8
EK 3813	Ø 3.8	13 mm	Ø 4.1	Ø 1.8
EK 3815	Ø 3.8	15 mm	Ø 4.1	Ø 1.8



REGULAR TAPERED SHAPE



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EK 4608	Ø 4.6	8.5 mm	Ø 4.1	Ø 2.6
EK 4610	Ø 4.6	10 mm	Ø 4.1	Ø 2.6
EK 4611	Ø 4.6	11.5 mm	Ø 4.1	Ø 2.6
EK 4613	Ø 4.6	13 mm	Ø 4.1	Ø 2.6
EK 4615	Ø 4.6	15 mm	Ø 4.1	Ø 2.6

