



# EXTERNAL HEXAGON

## SPEEDY MULTI PLATFORM


ES  
TAPERED



3.3  
NARROW



4.1  
REGULAR



5.0  
WIDE




IMPIANTO  
SENZA MOUNTER  
MOUNTERLESS


# ES 33

**S-TAPERED**  
external hexagon - full treatment

## NARROW SPEEDY SHAPE





CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 3310	Ø 3.3	10 mm	Ø 3.5	Ø 2.0
ES 3311	Ø 3.3	11,5 mm	Ø 3.5	Ø 2.0
ES 3313	Ø 3.3	13 mm	Ø 3.5	Ø 2.0
ES 3315	Ø 3.3	15 mm	Ø 3.5	Ø 2.0




**S-TAPERED HT**  
external hexagon - half treatment

**ES 33 HT**

NARROW SPEEDY SHAPE



CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 3310 HT	Ø 3.3	10 mm	Ø 3.5	Ø 2.0
ES 3311 HT	Ø 3.3	11,5 mm	Ø 3.5	Ø 2.0
ES 3313 HT	Ø 3.3	13 mm	Ø 3.5	Ø 2.0



# narrow platform prosthetics

20Ncm 



**Pilastro di Guarigione** - Titanio Gr. 5 - Anodizzato

**Healing Abutment** - Titanium Gr. 5 - Anodized

**SP 3HA3** H 3 mm - Narrow Gingival Profile

**SP 3HA5** H 5 mm - Narrow Gingival Profile

**DMX 1/2**

[Toronto Abutment Torx Driver](#)

**Analogo per Laboratorio** - Titanio Gr. 5 - Anodizzato

**Laboratory Analog** - Titanium Gr. 5 - Anodized



**SP 3AL** H 12 mm - Narrow Shape

10Ncm 



**Moncone Calcinabile** - Acrilico

**Castable Abutment** - Acrylic

**SP 3CA** Antirotaion Hex - Narrow Gingival Profile

**SP 3CR** Rotating - Narrow Gingival Profile

**SP VIT33** Vite di Ricambio - Spare Screw

**DMX 1/2**

[Toronto Abutment Torx Driver](#)

32Ncm 



**Moncone Titanio** - Titanio Gr. 5 - PVD TiN

**Titanium Abutment** - Titanio Gr. 5 - PVD TiN

**SP 3ST** 0° - Straight - Narrow Gingival Profile

**SP VIT33** Vite di Ricambio - Spare Screw

**DMX 1/2**

[Toronto Abutment Torx Driver](#)

32Ncm 



**Moncone Titanio** - Titanio Gr. 5 - PVD TiN

**Titanium Abutment** - Titanio Gr. 5 - PVD TiN

**SP 3AT15** 15° - Angled - Narrow Gingival Profile

**SP 3AT25** 25° - Angled - Narrow Gingival Profile

**SP VIT33** Vite di Ricambio - Spare Screw

**DMX 1/2**

[Toronto Abutment Torx Driver](#)

# narrow platform prosthetics

**Pilastro Temporaneo** - Titanio Gr. 5


**Temporary Abutment** - Titanium Gr. 5

**SP 3TA** Antirotaion Hex - Narrow Gingival Profile

**SP 3TR** Rotating - Narrow Gingival Profile

**SP VIT33** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver



**Transfer per Impronta Pick Up** - Titanio Gr. 5

**Pick Up Impression Coping** - Titanium Gr. 5

**SP 3IA** H 12mm - Narrow Gingival Profile

**SP VITL33** Vite di Ricambio Lunga - Long Spare Screw


**DMX 1/2** Toronto Abutment Torx Driver




# ES 38-40

**S-TAPERED**  
external hexagon - full treatment

## REGULAR SPEEDY SHAPE





CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 3808	Ø 3.8	8,5 mm	Ø 4.1	Ø 2.2
ES 3810	Ø 3.8	10 mm	Ø 4.1	Ø 2.2
ES 3811	Ø 3.8	11.5 mm	Ø 4.1	Ø 2.2
ES 3813	Ø 3.8	13 mm	Ø 4.1	Ø 2.2
ES 3815	Ø 3.8	15 mm	Ø 4.1	Ø 2.2
ES 4007	Ø 4.0	7 mm	Ø 4.1	Ø 2.5
ES 4008	Ø 4.0	8,5 mm	Ø 4.1	Ø 2.5
ES 4010	Ø 4.0	10 mm	Ø 4.1	Ø 2.5
ES 4011	Ø 4.0	11.5 mm	Ø 4.1	Ø 2.5
ES 4013	Ø 4.0	13 mm	Ø 4.1	Ø 2.5
ES 4015	Ø 4.0	15 mm	Ø 4.1	Ø 2.5
ES 4018	Ø 4.0	18 mm	Ø 4.1	Ø 2.5



**S-TAPERED HT**  
external hexagon - half treatment

**ES 40 HT**

REGULAR SPEEDY SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 4008	Ø 4.0	8,5 mm	Ø 4.1	Ø 2.5
ES 4010	Ø 4.0	10 mm	Ø 4.1	Ø 2.5
ES 4011	Ø 4.0	11.5 mm	Ø 4.1	Ø 2.5
ES 4013	Ø 4.0	13 mm	Ø 4.1	Ø 2.5
ES 4015	Ø 4.0	15 mm	Ø 4.1	Ø 2.5
ES 4018	Ø 4.0	18 mm	Ø 4.1	Ø 2.5

SOFT AND STD BONE - D2/D3/D4



HARD BONE - D1



# regular platform prosthetics

20Ncm 



**Pilastro di Guarigione** - Titanio Gr. 5 - Anodizzato

**Healing Abutment** - Titanium Gr. 5 - Anodized

**SP 4HA3** H 3 mm - Standard Gingival Profile

**SP 4HA5** H 5 mm - Standard Gingival Profile

**DMX 1/2**

Toronto Abutment Torx Driver



**Analogo per Laboratorio** - Titanio Gr. 5 - Anodizzato

**Laboratory Analog** - Titanium Gr. 5 - Anodized

**SP 4AL** H 12 mm - Standard Shape

10Ncm 



**Moncone Calcinabile Antirotazione** - Acrilico

**Antirootation Castable Abutment** - Acrylic

**SP 4CA** Antirootation Hex - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

10Ncm 



**Moncone Calcinabile Rotante** - Acrilico

**Rotating Castable Abutment** - Acrylic

**SP 4CR** Rotating - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

32Ncm 



**Moncone Titanio** - Titanio Gr. 5 - PVD TiN

**Titanium Abutment** - Titanio Gr. 5 - PVD TiN

**SP 4ST** 0° - Straight - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

# regular platform prosthetics

**Moncone Titanio** - Titanio Gr. 5 - PVD TiN

**Titanium Abutment** - Titanio Gr. 5 - PVD TiN

**SP 4AT15** 15° - Angled - Standard Gingival Profile

**SP 4AT25** 25° - Angled - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver



**Link per Incollaggio TiB** - Titanio Gr. 5

**Bonding Technique Link TiB** - Titanium Gr. 5

**EP 4STB** Link in Titanio Antirotazione - Compatibile TiB

**EP VIT** Vite di Ricambio - Spare Screw

**DM 1/2** Standard Hexagonal Driver



**Pilastro Temporaneo Antirotazione** - Titanio Gr. 5

**Antirotation Temporary Abutment** - Titanium Gr. 5

**SP 4TA** Antirotation Hex - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver



**Pilastro Temporaneo Rotante** - Titanio Gr. 5

**Rotating Temporary Abutment** - Titanium Gr. 5

**SP 4TR** Rotating - Standard Gingival Profile

**SP VIT40** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver




**Transfer per Impronta Pick Up** - Titanio Gr. 5

**Pick Up Impression Coping** - Titanium Gr. 5


**SP 4IA** H 12mm - Standard Gingival Profile

**SP VITL40** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver



# toronto universal prosthetics



## Pilastro Toronto 4 Dritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

- EP MS2      H 2 mm - Standard Gingival Profile
- EP MS3      H 3 mm - Standard Gingival Profile

**DCDM**      Toronto Straight Abutment Driver



## Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

- EP MA172    17° - H 2 mm - Standard Gingival Profile
- EP MA173    17° - H 3 mm - Standard Gingival Profile
- EP MA304    30° - H 4 mm - Standard Gingival Profile

EP VITMA    Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2**    Toronto Abutment Torx Driver



## Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

- EP MT      Rotating - Bonding Technique - Incollaggio Passivo
- EP MT-S    Rotating - Welding Technique - Saldatura Intraorale
- EP MV      Vite di Ricambio Toronto - Toronto Spare Screw
- DMX 1/2**    Toronto Abutment Torx Driver



# toronto universal prosthetics

10Ncm



**Pilastro Calcinabile Toronto** - Acrilico

**Toronto Castable Abutment** - Acrylic

**EP MC** Rotating - Standard Gingival Profile

**EP MV** Vite di Ricambio Toronto - Toronto Spare Screw

**DMX 1/2** [Toronto Abutment Torx Driver](#)

15Ncm



**Pilastro di Guarigione Toronto** - PE

**Toronto Healing Abutment** - PE

**EP MHA** H 4 mm - Standard Gingival Profile

**EP MHA6** H 6 mm - Standard Gingival Profile

**DMX 1/2** [Toronto Abutment Torx Driver](#)



15Ncm



**Transfer Toronto** - Titanio Gr. 5

**Toronto Impression Coping** - Titanium Gr. 5

**EP MI** Rotating - Standard Gingival Profile

**DMX 1/2** [Toronto Abutment Torx Driver](#)



**Analogo Toronto** - Titanio Gr. 5

**Toronto Analog** - Titanium Gr. 5

**EP MAL** H 10 mm - Universal Shape

15Ncm



15Ncm



**Viti di Ricambio Toronto** - Titanio Gr. 5

**Toronto Spare Screws** - Titanium Gr. 5

**EP MV** Prosthetic Spare Screw

**EP MVL** Impression Spare Screw

**DMX 1/2** [Toronto Abutment Torx Driver](#)


# ES 50

**S-TAPERED**  
external hexagon - full treatment


## WIDE SPEEDY SHAPE



- 0.5 mm Machined
- Tapered 1° Angled
- Screw Pitch 1.5 mm DT
- Tapered 3.5 mm




CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 5008	Ø 5.0	8,5 mm	Ø 5.0	Ø 3.5
ES 5010	Ø 5.0	10 mm	Ø 5.0	Ø 3.5
ES 5011	Ø 5.0	11.5 mm	Ø 5.0	Ø 3.5
ES 5013	Ø 5.0	13 mm	Ø 5.0	Ø 3.5



**S-TAPERED HT**  
external hexagon - half treatment

**ES 50 HT**

WIDE SPEEDY SHAPE



3.5 mm  
Machined

Tapered  
1° Angled

Screw Pitch  
1.6 mm DT


Tapered  
3.5 mm



CODE	IMPLANT	LENGTH	PLATFORM	APEX
ES 5008 HT	Ø 5.0	8,5 mm	Ø 5.0	Ø 3.5
ES 5010 HT	Ø 5.0	10 mm	Ø 5.0	Ø 3.5
ES 5011 HT	Ø 5.0	11.5 mm	Ø 5.0	Ø 3.5
ES 5013 HT	Ø 5.0	13 mm	Ø 5.0	Ø 3.5

**SOFT AND STD BONE - D2/D3/D4**

**HARD BONE - D1**



# wide platform prosthetics

20Ncm 



**Pilastro di Guarigione** - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

**SP 5HA3** H 3 mm - Large Gingival Profile

**SP 5HA5** H 5 mm - Large Gingival Profile

**DMX 1/2**

Toronto Abutment Torx Driver



**Analogo per Laboratorio** - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

**SP 5AL** H 12 mm - Large Shape

10Ncm 



**Moncone Calcinabile** - Acrilico

Castable Abutment - Acrylic

**SP 5CA** Antirotação Hex - Narrow Gingival Profile

**SP 5CR** Rotating - Large Gingival Profile

**SP VIT50** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

32Ncm 



**Moncone Titanio** - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

**SP 5ST** 0° - Straight - Large Gingival Profile

**SP VIT50** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

32Ncm 



**Pilastro Temporaneo con Esagono** - Titanio Gr. 5

Antirotation Hex Temporary Abutment - Titanium Gr. 5

**SP 5TA** Antirotação Hex - Narrow Gingival Profile

**SP VIT50** Vite di Ricambio - Spare Screw

**DMX 1/2**

Toronto Abutment Torx Driver

# wide platform prosthetics

**Pilastro Temporaneo senza Esagono - Titanio Gr. 5**

**Rotating Temporary Abutment - Titanium Gr. 5**

**SP 5TR** Rotating - Narrow Gingival Profile

**SP VIT50** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver



**Transfer per Impronta Pick Up - Titanio Gr. 5**

**Pick Up Impression Coping - Titanium Gr. 5**

**SP 5IA** H 12mm - Narrow Gingival Profile

**SP VITL50** Vite di Ricambio - Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver

