

narrow platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

SP 3HA3 H 3 mm - Narrow Gingival Profile

SP 3HA5 H 5 mm - Narrow Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

SP 3AL H 12 mm - Narrow Shape

10Ncm 



Moncone Calcinabile - Acrilico

Castable Abutment - Acrylic

SP 3CA Antirotation Hex - Narrow Gingival Profile

SP 3CR Rotating - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

SP 3ST 0° - Straight - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

SP 3AT15 15° - Angled - Narrow Gingival Profile

SP 3AT25 25° - Angled - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

narrow platform prosthetics

Pilastro Temporaneo - Titanio Gr. 5

Temporary Abutment - Titanium Gr. 5

SP 3TA Antirotation Hex - Narrow Gingival Profile

SP 3TR Rotating - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

SP 3IA H 12mm - Narrow Gingival Profile

SP VITL33 Vite di Ricambio Lunga - Long Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



regular platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

SP 4HA3 H 3 mm - Standard Gingival Profile

SP 4HA5 H 5 mm - Standard Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

SP 4AL H 12 mm - Standard Shape

10Ncm 



Moncone Calcinabile Antirotazione - Acrilico

Antirotation Castable Abutment - Acrylic

SP 4CA Antirotation Hex - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

10Ncm 



Moncone Calcinabile Rotante - Acrilico

Rotating Castable Abutment - Acrylic

SP 4CR Rotating - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

SP 4ST 0° - Straight - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

regular platform prosthetics

Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanio Gr. 5 - PVD TiN

SP 4AT15 15° - Angled - Standard Gingival Profile

SP 4AT25 25° - Angled - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



Link per Incollaggio TiB - Titanio Gr. 5

Bonding Technique Link TiB - Titanium Gr. 5

EP 4STB Link in Titanio Antirotazione - Compatibile TiB

EP VIT Vite di Ricambio - Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



Pilastro Temporaneo Antirotazione - Titanio Gr. 5

Antirotation Temporary Abutment - Titanium Gr. 5

SP 4TA Antirotation Hex - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



Pilastro Temporaneo Rotante - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

SP 4TR Rotating - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5


SP 4IA H 12mm - Standard Gingival Profile

SP VITL40 Vite di Ricambio - Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



toronto universal prosthetics



Pilastro Toronto 4 Diritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

EP MS2 H 2 mm - Standard Gingival Profile

EP MS3 H 3 mm - Standard Gingival Profile

DCDM

Toronto Straight Abutment Driver

 32Ncm



Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

EP MA172 17° - H 2 mm - Standard Gingival Profile

EP MA173 17° - H 3 mm - Standard Gingival Profile

EP MA304 30° - H 4 mm - Standard Gingival Profile

EP VITMA Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 **Toronto Abutment Torx Driver**

 32Ncm



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

EP MT Rotating - Bonding Technique - Incollaggio Passivo

EP MT-S Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 **Toronto Abutment Torx Driver**

15Ncm



toronto universal prosthetics

10Ncm



Pilastro Calcinabile Toronto - Acrilico

Toronto Castable Abutment - Acrylic

EP MC

Rotating - Standard Gingival Profile

EP MV

Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

15Ncm



Pilastro di Guarigione Toronto - PE

Toronto Healing Abutment - PE

EP MHA

H 4 mm - Standard Gingival Profile

EP MHA6

H 6 mm - Standard Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)



15Ncm



Transfer Toronto - Titanio Gr. 5

Toronto Impression Coping - Titanium Gr. 5

EP MI

Rotating - Standard Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)



Analogo Toronto - Titanio Gr. 5

Toronto Analog - Titanium Gr. 5

EP MAL

H 10 mm - Universal Shape

15Ncm



15Ncm



Viti di Ricambio Toronto - Titanio Gr. 5

Toronto Spare Screws - Titanium Gr. 5

EP MV

Prosthetic Spare Screw

EP MVL

Impression Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

wide platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

SP 5HA3 H 3 mm - Large Gingival Profile

SP 5HA5 H 5 mm - Large Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

SP 5AL H 12 mm - Large Shape

10Ncm 



Moncone Calcinabile - Acrilico

Castable Abutment - Acrylic

SP 5CA Antirotation Hex - Narrow Gingival Profile

SP 5CR Rotating - Large Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

SP 5ST 0° - Straight - Large Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Pilastro Temporaneo con Esagono - Titanio Gr. 5

Antirotation Hex Temporary Abutment - Titanium Gr. 5

SP 5TA Antirotation Hex - Narrow Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

wide platform prosthetics

Pilastro Temporaneo senza Esagono - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

SP 5TR Rotating - Narrow Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

SP 5IA H 12mm - Narrow Gingival Profile

SP VITL50 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

