



# EXTERNAL HEXAGON

## 4.1 REGULAR CONNECTION

SHORT  
CYLINDRICAL  
SHAPE



EC



IMPIANTO  
**CON MOUNTER**  
WITH MOUNTER


# EC 42

**CYLINDRIC**  
external hexagon - full treatment

## REGULAR CYLINDRICAL SHAPE



- 0.5 mm Machined
- Cylindric
- Screw Pitch 1.2 mm DT
- Tapered 3 mm




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 4206	Ø 4.2	6,0 mm	Ø 4.1	Ø 3.1

**SOFT AND STD BONE - D2/D3/D4**

**HARD BONE - D1**


WIDE CYLINDRICAL SHAPE



- 0.7 mm Machined
- Cylindric
- Screw Pitch 1.6 mm DT
- Tapered 3 mm




CODE	IMPLANT	LENGTH	PLATFORM	APEX
<b>EC 5006</b>	Ø 5.0	6.0 mm	Switch Ø 4.1	Ø 3.9




# EC 34 HT

**CYLINDRIC HT**  
external hexagon - half treatment

## NARROW CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3408 HT	Ø 3.4	8.5 mm	Ø 4.1	Ø 2.5
EC 3410 HT	Ø 3.4	10 mm	Ø 4.1	Ø 2.5
EC 3411 HT	Ø 3.4	11.5 mm	Ø 4.1	Ø 2.5
EC 3413 HT	Ø 3.4	13 mm	Ø 4.1	Ø 2.5



**CYLINDRIC HT**  
external hexagon - half treatment

**EC 38 HT**

REGULAR CYLINDRICAL SHAPE



3.5 mm  
Machined

Cylindric


Screw Pitch  
1.2 mm DT

Tapered  
3 mm




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 3808 HT	Ø 3.8	8.5 mm	Ø 4.1	Ø 2.7
EC 3810 HT	Ø 3.8	10 mm	Ø 4.1	Ø 2.7
EC 3811 HT	Ø 3.8	11.5 mm	Ø 4.1	Ø 2.7
EC 3813 HT	Ø 3.8	13 mm	Ø 4.1	Ø 2.7
EC 3815 HT	Ø 3.8	15 mm	Ø 4.1	Ø 2.7

SOFT AND STD BONE - D2/D3/D4




HARD BONE - D1




# EC 42 HT

**CYLINDRIC HT**  
external hexagon - half treatment

## REGULAR CYLINDRICAL SHAPE




CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 4208 HT	Ø 4.2	8.5 mm	Ø 4.1	Ø 3.1
EC 4210 HT	Ø 4.2	10 mm	Ø 4.1	Ø 3.1
EC 4211 HT	Ø 4.2	11.5 mm	Ø 4.1	Ø 3.1
EC 4213 HT	Ø 4.2	13 mm	Ø 4.1	Ø 3.1



**CYLINDRIC HT**  
external hexagon - half treatment

**EC 50 HT**

WIDE CYLINDRICAL SHAPE



3.5 mm  
Machined

Cylindric

Screw Pitch  
1.4 mm DT

Tapered  
3 mm



CODE	IMPLANT	LENGTH	PLATFORM	APEX
EC 5008 HT	Ø 5.0	8.5 mm	Switch Ø 4.1	Ø 3.9
EC 5010 HT	Ø 5.0	10 mm	Switch Ø 4.1	Ø 3.9
EC 5011 HT	Ø 5.0	11.5 mm	Switch Ø 4.1	Ø 3.9
EC 5013 HT	Ø 5.0	13 mm	Switch Ø 4.1	Ø 3.9

SOFT AND STD BONE - D2/D3/D4

HARD BONE - D1

