

EXTERNAL HEXAGON

4.1 REGULAR CONNECTION

UNIVERSAL PROSTHETICS



4.1 universal prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

EP 4HA2 H 2 mm - Standard Gingival Profile

EP 4HA4 H 4 mm - Standard Gingival Profile

EP 4HA6 H 6 mm - Standard Gingival Profile

DM 1/2 [Standard Hexagonal Driver](#)

20Ncm 



Pilastro di Guarigione Allargato - Titanio Gr. 5 - Anodizzato

Large Healing Abutment - Titanium Gr. 5 - Anodized

EP 5HA3 H 3 mm - **Large Gingival Profile**

EP 5HA5 H 5 mm - **Large Gingival Profile**

DM 1/2 [Standard Hexagonal Driver](#)




Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

EP 4AL H 12 mm - Universal Shape

10Ncm 



Moncone Calcinabile - Acrilico

Castable Abutment - Acrylic

EP 4CA Antirotation Hex - Standard Gingival Profile

EP 4CR Rotating - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

EP 4ST2 0° - Straight - Standard Gingival Profile

EP 5ST3 0° - Straight - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)

4.1 universal prosthetics

Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanio Gr. 5 - PVD TiN

EP 4AT15 15° - Angled - Standard Gingival Profile

EP 4AT25 25° - Angled - Standard Gingival Profile

EP 5AT15 15° - Angled - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Moncone UCLA Oro - Base Oro/Platino/Iridio

Ucla Gold Abutment - Gold/Platinum/Iridium Base

EP 4GA Antirotation Hex - Standard Gingival Profile

EP 5GA Antirotation Hex - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Moncone Provvisorio Estetico - PEEK Medicaale

Aesthetic Temporary Abutment - Medical PEEK

EP 4PA Antirotation Hex - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 25Ncm



Pilastro Temporaneo Antirotazione - Titanio Gr. 5

Antirotation Temporary Abutment - Titanium Gr. 5

EP 4TA Antirotation Hex - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Pilastro Temporaneo Rotante - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

EP 4TR Rotating - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



4.1 universal prosthetics

32Ncm 



Link per Incollaggio TiB - Titanio Gr. 5

Bonding Technique Link TiB - Titanium Gr. 5

EP 4STB

Link in Titanio Antirotazione - Compatibile TiB

EP VIT

Vite di Ricambio - Spare Screw

DM 1/2

Standard Hexagonal Driver

32Ncm 



Moncone Estetico - Zirconia

Esthetic Abutment - Zirconium

EP 4SZ

0° - Straight - Standard Gingival Profile

EP 4AZ15

15° - Angled - Standard Gingival Profile

EP VIT

Vite di Ricambio - Spare Screw

DM 1/2

Standard Hexagonal Driver

32Ncm 



DM Standard Hexagonal Driver

Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Abutment - Titanium Gr. 5 - PVD TiN

EP 4BA1

H 1 mm - Standard Gingival Profile - NORMO

EP 4BA2

H 2 mm - Standard Gingival Profile - NORMO

EP 4BA3

H 3 mm - Standard Gingival Profile - NORMO

EP 4BA4

H 4 mm - Standard Gingival Profile - NORMO

EP 4BA6

H 6 mm - Standard Gingival Profile - NORMO

32Ncm 



DMEQ Squared Equator Driver

Pilastro per Overdenture Equator - Titanio Gr. 5 - PVD TiN

Equator Overdenture Abutment - Titanium Gr. 5 - PVD TiN

EP 4EQ2

H 2 mm - Narrow Gingival Profile

EP 4EQ3

H 3 mm - Narrow Gingival Profile

EP 4EQ4


H 4 mm - Narrow Gingival Profile

EP 4EQ5

H 5 mm - Narrow Gingival Profile

EP 4EQ6

H 6 mm - Narrow Gingival Profile



 25Ncm

Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

EP X40

H 12mm - Standard Gingival Profile


EP VITL

Vite di Ricambio - Spare Screw

DM 1/2

Standard Hexagonal Driver

toronto 4.1 universal prosthetics



Pilastro Toronto 4 Dritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

- EP MS1** H 1 mm - Standard Gingival Profile
- EP MS2** H 2 mm - Standard Gingival Profile
- EP MS3** H 3 mm - Standard Gingival Profile

DCDM [Toronto Straight Abutment Driver](#)



32Ncm



Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

- EP MA172** 17° - H 2 mm - Standard Gingival Profile
- EP MA173** 17° - H 3 mm - Standard Gingival Profile
- EP MA304** 30° - H 4 mm - Standard Gingival Profile

EP VITMA Vite di Ricambio Toronto - Toronto Spare Screw
DMX 1/2 [Toronto Abutment Torx Driver](#)



32Ncm



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

- EP MT** Rotating - Bonding Technique - Incollaggio Passivo
- EP MT-S** Rotating - Welding Technique - Saldatura Intraorale
- EP MV** Vite di Ricambio Toronto - Toronto Spare Screw
- DMX 1/2** [Toronto Abutment Torx Driver](#)

15Ncm



toronto 4.1 universal prosthetics

10Ncm



Pilastro Calcinabile Toronto - Acrilico

Toronto Castable Abutment - Acrylic

EP MC Rotating - Standard Gingival Profile

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



15Ncm



Transfer Toronto - Titanio Gr. 5

Toronto Impression Coping - Titanium Gr. 5

EP MI Rotating - Standard Gingival Profile

EP MVL Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



Analogo Toronto - Titanio Gr. 5

Toronto Analog - Titanium Gr. 5

EP MAL H 10 mm - Universal Shape

15Ncm



Pilastro di Guarigione Toronto - PE

Toronto Healing Abutment - PE

EP MHA H 4 mm - Standard Gingival Profile

EP MHA6 H 6 mm - Standard Gingival Profile

DMX 1/2 [Toronto Abutment Torx Driver](#)

15Ncm



15Ncm



Viti di Ricambio Toronto - Titanio Gr. 5

Toronto Spare Screws - Titanium Gr. 5

EP MV Prosthetic Spare Screw

EP MVL Impression Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)