


CONICAL CONNECTION

NARROW PLATINUM

TAPERED
SHAPE



IK-P





IMPIANTO
SENZA MOUNTER
MOUNTERLESS


IK 34-P

NARROW K-PLATINUM
conical connection 12° - full treatment

TAPERED SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
IK 3408 P	Ø 3.4	8.5 mm	Ø 3.5	Ø 1.7
IK 3410 P	Ø 3.4	10 mm	Ø 3.5	Ø 1.7
IK 3411 P	Ø 3.4	11.5 mm	Ø 3.5	Ø 1.7
IK 3413 P	Ø 3.4	13 mm	Ø 3.5	Ø 1.7
IK 3415 P	Ø 3.4	15 mm	Ø 3.5	Ø 1.7




NARROW K-PLATINUM

conical connection 12° - full treatment

IK 38-P

TAPERED SHAPE



0.5 mm
Machined

Tapered 2°

Screw Pitch
0.8 mm

Tapered
6 mm



CODE	IMPLANT	LENGTH	NECK	APEX
IK 3808 P	Ø 3.8	8.5 mm	Ø 3.8	Ø 1.8
IK 3810 P	Ø 3.8	10 mm	Ø 3.8	Ø 1.8
IK 3811 P	Ø 3.8	11.5 mm	Ø 3.8	Ø 1.8
IK 3813 P	Ø 3.8	13 mm	Ø 3.8	Ø 1.8
IK 3815 P	Ø 3.8	15 mm	Ø 3.8	Ø 1.8

