


# CONICAL CONNECTION

## NARROW PLATINUM

### ACTIVE BLADE SHAPE



IA-P




IMPIANTO  
**SENZA MOUNTER**  
MOUNTERLESS


# IA 34-P

**NARROW A - PLATINUM**  
conical connection 12° - full treatment

## ACTIVE BLADE SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
IA 3408 P	Ø 3.4	8.5 mm	Ø 3.5	Ø 1.8
IA 3410 P	Ø 3.4	10 mm	Ø 3.5	Ø 1.8
IA 3411 P	Ø 3.4	11.5 mm	Ø 3.5	Ø 1.8
IA 3413 P	Ø 3.4	13 mm	Ø 3.5	Ø 1.8
IA 3415 P	Ø 3.4	15 mm	Ø 3.5	Ø 1.8




# NARROW A - PLATINUM

conical connection 12° - full treatment

# IA 38-P

## ACTIVE BLADE SHAPE



1.5 mm  
Machined

Cylindric


Screw Pitch  
1.2 mm

Tapered  
4 mm



CODE	IMPLANT	LENGTH	NECK	APEX
IA 3808 P	Ø 3.8	8.5 mm	Ø 3.8	Ø 2.1
IA 3810 P	Ø 3.8	10 mm	Ø 3.8	Ø 2.1
IA 3811 P	Ø 3.8	11.5 mm	Ø 3.8	Ø 2.1
IA 3813 P	Ø 3.8	13 mm	Ø 3.8	Ø 2.1
IA 3815 P	Ø 3.8	15 mm	Ø 3.8	Ø 2.1

### SOFT AND STD BONE - D2/D3/D4



### HARD BONE - D1

