



# CONICAL CONNECTION

## NARROW PLATINUM

SPEEDY  
SHAPE



ES-P




IMPIANTO  
**SENZA MOUNTER**  
MOUNTERLESS


# ES 34-P

**NARROW S-PLATINUM**  
conical connection 12° - full treatment

## SPEEDY SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
ES 3408 P	Ø 3.4	8.5 mm	Ø 3.5	Ø 1.8
ES 3410 P	Ø 3.4	10 mm	Ø 3.5	Ø 1.8
ES 3411 P	Ø 3.4	11.5 mm	Ø 3.5	Ø 1.8
ES 3413 P	Ø 3.4	13 mm	Ø 3.5	Ø 1.8
ES 3415 P	Ø 3.4	15 mm	Ø 3.5	Ø 1.8




# NARROW S - PLATINUM

conical connection 12° - full treatment

# ES 38-P

## ACTIVE BLADE SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
ES 3808 P	Ø 3.8	8.5 mm	Ø 3.8	Ø 2.0
ES 3810 P	Ø 3.8	10 mm	Ø 3.8	Ø 2.0
ES 3811 P	Ø 3.8	11.5 mm	Ø 3.8	Ø 2.0
ES 3813 P	Ø 3.8	13 mm	Ø 3.8	Ø 2.0
ES 3815 P	Ø 3.8	15 mm	Ø 3.8	Ø 2.0

