



# CONICAL CONNECTION

## REGULAR PLATINUM

TAPERED  
SHAPE



IK-P




IMPIANTO  
**SENZA MOUNTER**  
MOUNTERLESS


# IK 42-P

**REGULAR K - PLATINUM**  
conical connection 12° - full treatment

## TAPERED SHAPE





CODE	IMPLANT	LENGTH	NECK	APEX
IK 4207 P	Ø 4.2	7.0 mm	Ø 4.0	Ø 2.2
IK 4208 P	Ø 4.2	8.5 mm	Ø 4.0	Ø 2.2
IK 4210 P	Ø 4.2	10 mm	Ø 4.0	Ø 2.2
IK 4211 P	Ø 4.2	11.5 mm	Ø 4.0	Ø 2.2
IK 4213 P	Ø 4.2	13 mm	Ø 4.0	Ø 2.2
IK 4215 P	Ø 4.2	15 mm	Ø 4.0	Ø 2.2



**REGULAR K - PLATINUM**  
conical connection 12° - full treatment

**IK 50-P**

TAPERED SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
IK 5007 P	Ø 5.0	7.0 mm	Ø 4.8	Ø 3.3
IK 5008 P	Ø 5.0	8.5 mm	Ø 4.8	Ø 3.3
IK 5010 P	Ø 5.0	10 mm	Ø 4.8	Ø 3.3
IK 5011 P	Ø 5.0	11.5 mm	Ø 4.8	Ø 3.3
IK 5013 P	Ø 5.0	13 mm	Ø 4.8	Ø 3.3

