



CONICAL CONNECTION

REGULAR PLATINUM

ACTIVE BLADE SHAPE



IA-P




IMPIANTO
SENZA MOUNTER
MOUNTERLESS


IA 42-P

REGULAR A - PLATINUM
conical connection 12° - full treatment

ACTIVE BLADE SHAPE




CODE	IMPLANT	LENGTH	NECK	APEX
IA 4208 P	Ø 4.2	8.5 mm	Ø 4.0	Ø 2.2
IA 4210 P	Ø 4.2	10 mm	Ø 4.0	Ø 2.2
IA 4211 P	Ø 4.2	11.5 mm	Ø 4.0	Ø 2.2
IA 4213 P	Ø 4.2	13 mm	Ø 4.0	Ø 2.2
IA 4215 P	Ø 4.2	15 mm	Ø 4.0	Ø 2.2



REGULAR A - PLATINUM
conical connection 12° - full treatment

IA 50-P

ACTIVE BLADE SHAPE



CODE	IMPLANT	LENGTH	NECK	APEX
IA 5008 P	Ø 5.0	8.5 mm	Ø 4.8	Ø 3.0
IA 5010 P	Ø 5.0	10 mm	Ø 4.8	Ø 3.0
IA 5011 P	Ø 5.0	11.5 mm	Ø 4.8	Ø 3.0
IA 5013 P	Ø 5.0	13 mm	Ø 4.8	Ø 3.0

