



# CONICAL CONNECTION

## REGULAR PLATINUM

SPEEDY  
SHAPE



ES-P




IMPIANTO  
**SENZA MOUNTER**  
MOUNTERLESS


# ES 42-P

**REGULAR S - PLATINUM**  
conical connection 12° - full treatment

## SPEEDY SHAPE




CODE	IMPLANT	LENGTH	NECK	APEX
ES 4208 P	Ø 4.2	8.5 mm	Ø 4.0	Ø 2.2
ES 4210 P	Ø 4.2	10 mm	Ø 4.0	Ø 2.2
ES 4211 P	Ø 4.2	11.5 mm	Ø 4.0	Ø 2.2
ES 4213 P	Ø 4.2	13 mm	Ø 4.0	Ø 2.2
ES 4215 P	Ø 4.2	15 mm	Ø 4.0	Ø 2.2



**REGULAR S - PLATINUM**  
conical connection 12° - full treatment

**ES 50-P**

**SPEEDY SHAPE**



- 0.5 mm Machined
- Tapered 2° Angled
- Screw Pitch 1.6 mm DT
- Tapered 4 mm



CODE	IMPLANT	LENGTH	NECK	APEX
ES 5008 P	Ø 5.0	8.5 mm	Ø 4.8	Ø 3.5
ES 5010 P	Ø 5.0	10 mm	Ø 4.8	Ø 3.5
ES 5011 P	Ø 5.0	11.5 mm	Ø 4.8	Ø 3.5
ES 5013 P	Ø 5.0	13 mm	Ø 4.8	Ø 3.5

