




pterygoids and nasal implants



TAPERED SHAPE IK 42



CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 4218	Ø 4.2	18 mm	Ø 4.0	Ø 2.2
IK 4220	Ø 4.2	20 mm	Ø 4.0	Ø 2.2
IK 4222	Ø 4.2	22 mm	Ø 4.0	Ø 2.2
IK 4225	Ø 4.2	25 mm	Ø 4.0	Ø 2.2



pterygoids and nasal implants



pterygoids and nasal implants

