




Kit Chirurgical

REVOLUTION




REVOLUTION



REVOLUTION

59 STRUMENTI ROTANTI

MUCOTOMI – BONE MILL – FRESE CALIBRATE





FRESE CORTICALI




FRESA PIN + PINS



ACCESSORI




MONTATORI




MUCOTOMI



FRESE CHIRURGICHE






**SONDA MILLIMETRATA
BOTTONUTA**



**CRICCHETTO
DINAMOMETRICO**



Frese e Boccole

REVOLUTION

FRESE S

Ø4.15 mm



Small

FRESE L

Ø5.5 mm



Large




Small




Large



RAPPORTO IMPIANTO - BOCCOLA




Boccole SMALL $\varnothing 4.15$ mm



IMPIANTI : $\varnothing 3.0$ - $\varnothing 3.4$ - $\varnothing 3.8$



Boccole LARGE $\varnothing 5.5$ mm




IMPIANTI : $\varnothing 3.0$ - $\varnothing 3.4$ - $\varnothing 3.8$ - $\varnothing 4.2$ - $\varnothing 4.6$ - $\varnothing 5.0$



Mounter

REVOLUTION

MOUNTER **MTA** PREASSEMBLATO



MOUNTER CALIBRATO STOPPATO

ALTERNATIVA CHIRURGICA

Ø4.15



Ø5.5



● REPERE DI FASE

● STOP

● GUIDA

● SCARICO


● FACCIA ESAGONO
IN FASE

MOUNTER CALIBRATO




BASE IN TITANIO RIVESTITA DLC


RIMOZIONE DEL MOUNTER STANDARD



1



2



3