

CONICAL CONNECTION

REGULAR PLATINUM

REGULAR
PROSTHETICS



regular CCP prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

IP 4HA3-P H 3 mm - Conical Gingival Profile

IP 4HA4-P H 4 mm - Conical Gingival Profile

IP 4HA5-P H 5 mm - Conical Gingival Profile

DMX 1/2

Toronto Abutment Torx Driver


Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized



IP 4AL-P H 12 mm - Universal Shape

32Ncm 



Moncone UCLA Cobalto Cromo - CoCr

UCLA Abutment Cobalt Chrome - CoCr

IP 4CCA-P Antirrotation Hex - Conical Gingival Profile H1,5mm

IP 4CCR-P Rotating - Conical Gingival Profile H1,5mm

SP VIT40

Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

32Ncm 



Moncone Titanio Diritto - Titanio Gr. 5

Straight Titanium Abutment - Titanium Gr. 5

IP 4ST1-P 0° - Straight - Conical Gingival Profile H1.5mm

IP 4ST3-P 0° - Straight - Conical Gingival Profile H3.0mm

IP 4ST4-P 0° - Straight - Conical Gingival Profile H4.0mm

SP VIT40

Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

32Ncm 



Moncone Titanio Preangolato 15° - Titanio Gr. 5

15° Angled Titanium Abutment - Titanium Gr. 5

IP 4AT151-P 15° - Angled - Conical Gingival Profile H1.5mm

IP 4AT153-P 15° - Angled - Conical Gingival Profile H3.0mm

IP 4AT154-P 15° - Angled - Conical Gingival Profile H4.0mm

SP VIT40

Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

regular CCP prosthetics

Moncone Titanio Preangolato 25° - Titanio Gr. 5

25° Angled Titanium Abutment - Titanium Gr. 5


IP 4AT251-P 25° - Angled - Conical Gingival Profile H1.5mm

IP 4AT253-P 25° - Angled - Conical Gingival Profile H3.0mm

IP 4AT254-P 25° - Angled - Conical Gingival Profile H4.0mm

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo da Incollaggio Antirotazione - Ti Gr. 5

Antirotation Hex Temporary Glueing Abutment - Ti Gr. 5

IP 4TA1-P Antirotation Hex - Standard Gingival Profile H 1.5mm

IP 4TA3-P Antirotation Hex - Standard Gingival Profile H 3.0mm

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo da Incollaggio Rotante - Titanio Gr. 5

Rotating Temporary Glueing Abutment - Titanium Gr. 5

IP 4TR1-P Rotating - Conical Gingival Profile H 1.5mm

IP 4TR3-P Rotating - Conical Gingival Profile H 3.0mm

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Abutment - Titanium Gr. 5 - PVD TiN

IP 4BA1-P H 1 mm - Conical Gingival Profile - NORMO

IP 4BA2-P H 2 mm - Conical Gingival Profile - NORMO

IP 4BA3-P H 3 mm - Conical Gingival Profile - NORMO

IP 4BA4-P H 4 mm - Conical Gingival Profile - NORMO

DMX 1/2 Toronto Abutment Torx Driver

 32Ncm



Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

IP 4IA1-P H 12mm - Internal Hexagon H1.5mm


IP 4IA3-P H 12mm - Internal Hexagon H3.0mm

IP 4IA4-P H 12mm - Internal Hexagon H4.0mm


SP VITL40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

25Ncm 



toronto regular CCP prosthetics



32Ncm 



Pilastro Toronto Diritto - Titanio Gr. 5

Toronto Straight Abutment - Titanium Gr. 5

IP 4MS2-P H 2.0 mm - Standard Gingival Profile


IP 4MS3-P H 3.0 mm - Standard Gingival Profile

IP 4MS4-P H 4.0 mm - Standard Gingival Profile

DCDM

[Toronto Straight Abutment Driver](#)

32Ncm 



Pilastro Toronto Angolato 17° - Titanio Gr. 5

Toronto 17° Angled Abutment - Titanium Gr. 5

IP 4MA172-P 17° - H 2 mm - Standard Gingival Profile

IP 4MA173-P 17° - H 3 mm - Standard Gingival Profile


EP VITMA

Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2

[Toronto Abutment Torx Driver](#)

32Ncm 



Pilastro Toronto Angolato 30° - Titanio Gr. 5

Toronto 30° Angled Abutment - Titanium Gr. 5

IP 4MA303-P 30° - H 3 mm - Standard Gingival Profile

IP 4MA304-P 30° - H 4 mm - Standard Gingival Profile

DMX 1/2

[Toronto Abutment Torx Driver](#)

toronto regular CCP prosthetics

Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

EP MT Rotating - Bonding Technique - Incollaggio Passivo

EP MT-S Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Calcinabile Toronto - Acrilico

Toronto Castable Abutment - Acrylic

EP MC Rotating - Standard Gingival Profile

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Componenti Toronto - Per Impronta e Modello in Gesso

Toronto Components - For Impression and Gypsum Model

EP MI Toronto Impression Abutment - Transfer Impronta

EP MAL Toronto Analog - Analogo da Laboratorio

EP MVL Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro di Guarigione Toronto - PE

Toronto Healing Abutment - PE

EP MHA H 4 mm - Standard Gingival Profile

EP MHA6 H 6 mm - Standard Gingival Profile

DMX 1/2 Toronto Abutment Torx Driver



Viti di Ricambio Toronto - Titanio Gr. 5

Toronto Spare Screws - Titanium Gr. 5

EP MV Prosthetic Spare Screw

EP MVL Impression Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

