

OVERDENTURE & PROSTHETICS

ATTACHMENTS





Kit OT Equator Esagono Interno - Pilastro + Contenitore + Cappette

Kit OT Equator Internal Hexagon - Abutment + Box + Caps

IP3EQ2	Kit OT Equator h2 mm Ø3.5	●
IP3EQ3	Kit OT Equator h3 mm Ø3.5	●
IP3EQ4	Kit OT Equator h4 mm Ø3.5	●
IP3EQ5	Kit OT Equator h5 mm Ø3.5	●
IP3EQ6	Kit OT Equator h6 mm Ø3.5	●
IP4EQ1	Kit OT Equator h1 mm Ø4.0 CPS	●
IP4EQ2	Kit OT Equator h2 mm Ø4.0 CPS	●
IP4EQ3	Kit OT Equator h3 mm Ø4.0 CPS	●
IP4EQ4	Kit OT Equator h4 mm Ø4.0 CPS	●
IP4EQ5	Kit OT Equator h5 mm Ø4.0 CPS	●
IP4EQ6	Kit OT Equator h6 mm Ø4.0 CPS	●



3.3

Kit OT Equator Esagono Esterno - Pilastro + Contenitore + Cappette

Kit OT Equator External Hexagon - Abutment + Box + Caps

EP3EQ2	Kit OT Equator h2 mm Ø3.3	●
EP3EQ3	Kit OT Equator h3 mm Ø3.3	●
EP3EQ4	Kit OT Equator h4 mm Ø3.3	●
EP3EQ5	Kit OT Equator h5 mm Ø3.3	●
EP3EQ6	Kit OT Equator h6 mm Ø3.3	●



4.1

Kit OT Equator Esagono Esterno - Pilastro + Contenitore + Cappette

Kit OT Equator External Hexagon - Abutment + Box + Caps

EP4EQ2	Kit OT Equator h2 mm Ø4.1	●●
EP4EQ3	Kit OT Equator h3 mm Ø4.1	●●
EP4EQ4	Kit OT Equator h4 mm Ø4.1	●●
EP4EQ5	Kit OT Equator h5 mm Ø4.1	●●
EP4EQ6	Kit OT Equator h6 mm Ø4.1	●●



overdenture

Kit OT Equator Esagono Esterno - Pilastro + Contenitore + Cappette

Kit OT Equator External Hexagon - Abutment + Box + Caps

EP5EQ2 Kit OT Equator h2 mm Ø5.0 ●

EP5EQ3 Kit OT Equator h3 mm Ø5.0 ●

EP5EQ4 Kit OT Equator h4 mm Ø5.0 ●

EP5EQ5 Kit OT Equator h5 mm Ø5.0 ●

EP5EQ6 Kit OT Equator h6 mm Ø5.0 ●

5.0



Kit OT Equator Connessione Conica - Pilastro + Contenitore + Cappette

Kit OT Equator Conical Connection - Abutment + Box + Caps

IP3EQ1-P Kit OT Equator h1 mm Conical C. - NP ●

IP3EQ2-P Kit OT Equator h2 mm Conical C. - NP ●

IP3EQ3-P Kit OT Equator h3 mm Conical C. - NP ●

IP3EQ4-P Kit OT Equator h4 mm Conical C. - NP ●

IP3EQ5-P Kit OT Equator h5 mm Conical C. - NP ●

IP3EQ6-P Kit OT Equator h6 mm Conical C. - NP ●

NP



Kit OT Equator Connessione Conica - Pilastro + Contenitore + Cappette

Kit OT Equator Conical Connection - Abutment + Box + Caps

IP4EQ1-P Kit OT Equator h1 mm Conical C. - RP ●

IP4EQ2-P Kit OT Equator h2 mm Conical C. - RP ●

IP4EQ3-P Kit OT Equator h3 mm Conical C. - RP ●


IP4EQ4-P Kit OT Equator h4 mm Conical C. - RP ●

IP4EQ5-P Kit OT Equator h5 mm Conical C. - RP ●

IP4EQ6-P Kit OT Equator h6 mm Conical C. - RP ●

RP





Femmine per Overdenture Normo - Ritenzioni Rhein

Normo Overdenture Caps - Rhein Caps

TF 10	Teflon Bianca 1300 gr.
TF 08	Teflon Rosa 800 gr.
TF 06	Teflon Gialla 600 gr.
TF 04	Teflon Verde 400 gr.
TF A	Sottodimensionata - Teflon Arancione 350 gr.
TF G	Sottodimensionata - Teflon Gialla 500 gr.
TF R	Sottodimensionata - Teflon Rosa 900 gr.
TF V	Sottodimensionata - Teflon Verde 1300 gr.
CO 01	Metal Box - Stainless Steel
CO 02	Metal Box - Titanium



Femmine per Overdenture Micro - Ritenzioni Rhein

Micro Overdenture Caps - Rhein Caps

TFM 10	Teflon Bianca 1100 gr.
TFM 08	Teflon Rosa 800 gr.
TFM 06	Teflon Gialla 450 gr.
TFM 04	Teflon Verde 200 gr.
TFM A	Sottodimensionata - Teflon Arancione 200 gr.
TFM V	Sottodimensionata - Teflon Verde 1100 gr.
COM 01	Metal Box - Stainless Steel
COM 02	Metal Box - Titanium



Femmine per Overdenture Equator - Ritenzioni Rhein


Equator Overdenture Caps - Rhein Caps

TF Q1	Teflon Viola Forte
TF Q2	Teflon Bianco Normale
TF Q3	Teflon Rosa Morbida
TF Q4	Teflon Giallo Extra
TF Q5	Teflon Nero da Laboratorio
SBOX EQ	Smart Box - Autoparallelizzante Cont. + 5 Teflon
COQ 01	Metal Box - Stainless Steel
COQ 02	Metal Box - Titanium




overdenture

AN D	Anelli Direzionali 0° - 7° - 14° - Conf 6pz.
F20 SN	Sfera Normo Filettata - Passo Filetto M2.0
F16 SN	Sfera Normo Filettata - Passo Filetto M1.6
F20 SM	Sfera Micro Filettata - Passo Filetto M2.0
F16 SM	Sfera Micro Filettata - Passo Filetto M1.6
SGU SN	Spaziatore Sfera Normo per Guaina Filettata
SGU SM	Spaziatore Sfera Micro per Guaina Filettata
GU 16	Guaina Filettata da Incollaggio - Passo Filetto M1.6



SC N	Sfera Calcinabile Normo Verde - 4pz.
SC M	Sfera Calcinabile Micro Rossa - 4pz.
DPI	Dischetto Protettivo per Incollaggio Trasparente - 10pz.
SC 25	Sfera Cava da Incollaggio Normo - Per sfere usurate - 2pz.
SC 18	Sfera Cava da Incollaggio Micro - Per sfere usurate - 2pz.
BA 01	Barra tipo Dolder Calcinabile
BA 02	Ritenzione in Teflon Rosa per Barra
BA 03	Ritenzione in Teflon Gialla per Barra



BCA N	Barra Calcinabile Angolare Normo Verde - 2pz.
BCD N	Barra Calcinabile Diritta Normo Verde - 2pz.
BCA M	Barra Calcinabile Angolare Micro Rossa - 2pz.
BCD M	Barra Calcinabile Diritta Micro Rossa - 2pz.
F20 EQ	Equator Filettato - Passo Filetto M2.0
F16 EQ	Equator Filettato - Passo Filetto M1.6
SGU EQ	Equator Spaziatore per Guaina Filettata
AL EQ	Analogo per Laboratorio Equator






Kit LOCATOR® Esagono Interno - Pilastro + Contenitore + Caps

Kit LOCATOR® Internal Hexagon - Abutment + Box + Caps

- IP4LC1 Kit LOCATOR h1 mm - Ø4.0 CPS ●
- IP4LC2 Kit LOCATOR h2 mm - Ø4.0 CPS ●
- IP4LC3 Kit LOCATOR h3 mm - Ø4.0 CPS ●
- IP4LC4 Kit LOCATOR h4 mm - Ø4.0 CPS ●
- IP4LC5 Kit LOCATOR h5 mm - Ø4.0 CPS ●
- IP4LC6 Kit LOCATOR h6 mm - Ø4.0 CPS ●

LOCATOR®

4.1



Kit LOCATOR® Esagono Esterno - Pilastro + Contenitore + Caps

Kit LOCATOR® External Hexagon - Abutment + Box + Caps

- EP4LC1 Kit LOCATOR h1 mm - Ø4.1 ●●
- EP4LC2 Kit LOCATOR h2 mm - Ø4.1 ●●
- EP4LC3 Kit LOCATOR h3 mm - Ø4.1 ●●
- EP4LC4 Kit LOCATOR h4 mm - Ø4.1 ●●
- EP4LC5 Kit LOCATOR h5 mm - Ø4.1 ●●

LOCATOR®

NP



Kit LOCATOR® Connessione Conica - Pilastro + Contenitore + Caps

Kit LOCATOR® Conical Connection - Abutment + Box + Caps

- IP3LC1-P Kit LOCATOR h1.5 mm - Conical C. - NP ●
- IP3LC2-P Kit LOCATOR h2.5 mm - Conical C. - NP ●
- IP3LC5-P Kit LOCATOR h5.0 mm - Conical C. - NP ●

LOCATOR®

overdenture

Kit LOCATOR® Connessione Conica - Pilastro + Contenitore + Cappette

Kit LOCATOR® Conical Connection - Abutment + Box + Caps

IP4LC1-P Kit LOCATOR h1.5 mm - Conical C. - RP ●

IP4LC2-P Kit LOCATOR h2.5 mm - Conical C. - RP ●

IP4LC5-P Kit LOCATOR h5.0 mm - Conical C. - RP ●



LOCATOR®

Accessori e Ricambi LOCATOR®

LOCATOR® Accessories and Spare Parts

W2-TSLC Transfer Ti + Cappetta Nera (per laboratorio)

W2-COLC1 Contenitore Ti + Cappetta Nera (per laboratorio)

W2-NBOXLC Kit Cappette + Contenitore Ti

W2-TFLC50 Cappetta Ritentiva Trasparente - 2268 g 4 pz.

W2-TFLC30 Cappetta Ritentiva Rosa - 1361 g 4 pz.

W2-TFLC15 Cappetta Ritentiva Blu - 680 g 4 pz.

W2-ALLC Analogo in Titanio

W2-TFLCD40 Cappetta Ritentiva Verde - 1814 g 4 pz.


W2-TFLCD20 Cappetta ritentiva Arancione - 907 g 4 pz.

W2-TFLCD15 Cappetta ritentiva Rossa - 680 g 4 pz.

W2-DMLC Avvitatore manuale da Cricchetto (attacco Straumann®)

W2-DCLC Avvitatore da Contrangolo

W2-IELC Inseritore / Estrattore Cappette



LOCATOR®