


INSTRUMENTS & TOOLS


IMPLANT TOOLS

DRIVERS



instruments and tools

CACCIAVITI MANUALI - MANUAL SCREW DRIVERS




DM 01 Hexagonal 1.27 - Short - 12mm



DM 02 Hexagonal 1.27 - Medium - 17mm

DM 03 Hexagonal 1.27 - Long - 27mm



DMX 01 Torx - Short - 15mm

DMX 02 Torx - Long - 20mm



DM T6DS T6 Angled Torx Driver - 15mm



DMEQ Equator Squared - Short - 15mm

CACCIAVITI DA CONTRANGOLO - CONTRA-ANGLE SCREW DRIVERS



DC 01 Hexagonal 1.27 - Short - 15mm



DC 02 Hexagonal 1.27 - Medium - 20mm

DC 03 Hexagonal 1.27 - Long - 30mm



DC T6DS T6 Angled Torx Driver - 20mm



DCEQ Equator Squared - Medium - 20mm

DCX 01 Torx - Short - 15mm



DCX 02 Torx - Long - 20mm

instruments and tools

CACCIAVITI DA CRICCHETTO - TORQUE CONTROLLER DRIVERS

DCD 01 Hexagonal 1.27 - Short - 15mm



DCD 02 Hexagonal 1.27 - Medium - 20mm

DCD 03 Hexagonal 1.27 - Long - 30mm



DCDX 01 Torx - Short - 15mm



DCDX 02 Torx - Long - 20mm


DCDM Hexagonal - Monconi Toronto Diritti



AVVITATORI MANUALI - FINGER MANUAL DEVICE

DGM Digitale Manuale - Manual Implant Driver

DKX Avvitatore Manuale per RA - RA Manual Driver



CRD-A Ghiera per CRD con attacco RA - CRD Ferrule with RA adapter

CRD-G Ghiera per CRD con attacco PC - CRD Ferrule with PC adapter



NB Controllare compatibilità vecchio e nuovo ed in alternativa catalogo REPLICA

AVVITATORE PER MINI IMPIANTI - MINI IMPLANT DRIVER

PCM Cricchetto - Ratchet **Mini Implant** Driver



ACM Contrangolo - Contrangle **Mini Implant** Driver



instruments and tools

PROLUNGHE DA CRICCHETTO - RATCHET IMPLANT DRIVERS



PC 01R Short - Hexagonal Implant Driver

I E

PC 02R Medium - Hexagonal Implant Driver



PC 03R Long - Hexagonal Implant Driver

PROLUNGHE DA CRICCHETTO - RATCHET IMPLANT DRIVERS



PC ES1 Medium - **Star Grip** Driver for **ES 40**

ES

PC ES2 Long - **Star Grip** Driver for **ES 40**



PC ES3S Short - **Star Grip** Driver for **ES 33**

PC ES3L Long - **Star Grip** Driver for **ES 33**



PC ES5S Short - **Star Grip** Driver for **ES 50**

PC ES5L Long - **Star Grip** Driver for **ES 50**

PROLUNGHE DA CRICCHETTO - RATCHET IMPLANT DRIVERS



PC P3S Short - Conical Grip Driver **Platinum**

P

PC P3M Medium - Conical Grip Driver **Platinum**

PC P3L Long - Conical Grip Driver **Platinum**

PC P4S Short - Conical Grip Driver **Platinum**



PC P4M Medium - Conical Grip Driver **Platinum**

PC P4L Long - Conical Grip Driver **Platinum**

instruments and tools

AVVITATORI DA CONTRANGOLO - CONTRA-ANGLE IMPLANT DRIVERS

AC 01 Short - Hexagonal Implant Driver



AC 02 Medium - Hexagonal Implant Driver



I E

AVVITATORI DA CONTRANGOLO - CONTRA-ANGLE IMPLANT DRIVERS

AC ES1 Short - Star Grip Driver for **ES 40**



AC ES2 Medium - Star Grip Driver for **ES 40**

AC ES301 Short - Star Grip Driver for **ES 33**



AC ES302 Medium - Star Grip Driver for **ES 33**

AC ES501 Short - Star Grip Driver for **ES 50**



AC ES502 Medium - Star Grip Driver for **ES 50**

ES

AVVITATORI DA CONTRANGOLO - CONTRA-ANGLE IMPLANT DRIVERS

AC P3S Short - Conical Grip Driver **Platinum**

AC P3M Medium - Conical Grip Driver **Platinum**



AC P3L Long - Conical Grip Driver **Platinum**

AC P4S Short - Conical Grip Driver **Platinum**

AC P4M Medium - Conical Grip Driver **Platinum**



AC P4L Long - Conical Grip Driver **Platinum**

P

instruments and tools

Perni di Parallelismo - Titanio Gr. 5



Paralleling Pins - Titanium Gr. 5

PIN Double diameter - Small 1.8mm / Large 3mm - 0°

PIN17 Double diameter - Small 1.8mm / Large 3mm - 17°

PIN32 Double diameter - Small 1.8mm / Large 3mm - 32°

PIN K Kit Misuratori di Parallelismo Calibrati (10 PIN + BOX)





Osseo Espansori e Compattatori - Titanio Gr. 5

Bone Expanders and Compactors - Titanium Gr. 5

EXP 4 Progressive diameters

KIT 005 Expander Surgical Kit



A 1.3	A 1.6	A 1.9	A 2.4
B 2.1	B 2.8	B 3.5	B 4.3
C 2.4	C 3.1	C 3.8	C 4.6

Stop per Frese Resista - Titanio Gr. 5

Resista Drill Stoppers - Titanium Gr. 5

KIT ST 6 Titanium Drill Stoppers (6.0/7.0/8.5/10/11.5/13)

STP KIT 11 Titanium Drill Stoppers (1.0/2.0/3.0/4.0/5.0/6.0/7.0/8.5/10/11.5/13)


ST Titanium Drill Stopper 1 pz. (Es. ST1, ST2, ST3)



instruments and tools

Avvitatori Manuali - Acciaio Inox

Manual Driver - Stainless Steel




CRND Cricchetto non Dinamometrico - No Torque Control



CRD2 Cricchetto Dinamometrico Chirurgico 0/70 Ncm - Torque Control


CRDP Cricchetto Dinamometrico Protesico 0/35 Ncm - Torque Control




CDL Avvitatore Diritto Manuale - Long Manual Driver

Sonde e Misuratori - Acciaio Inox

Probes and Meters - Stainless Steel



ST PRO SL1 Sonda Compattatore Doppia - Sinus Lift System



SOM Sonda Chirurgica in Titanio - Titanium Surgical Probe



RIG Righello Chirurgico in Titanio - Titanium Surgical Ruler