


INTERNAL HEXAGON


WIDE PLATFORM

DOUBLE CONNECTION

X-LARGE
TAPERED SHAPE



IK XL




IMPIANTO
SENZA MOUNTER
MOUNTERLESS

IK 60

K-TAPERED
internal hexagon - full treatment

X-LARGE TAPERED SHAPE



CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 6007	Ø 6.0	7.0 mm	Ø 5.0	Ø 4.5
IK 6008	Ø 6.0	8.5 mm	Ø 5.0	Ø 4.5
IK 6010	Ø 6.0	10 mm	Ø 5.0	Ø 4.5


ONLY FOR SOFT AND STD BONE - D3/D4



K-TAPERED
internal hexagon - full treatment

IK 80

X-LARGE TAPERED SHAPE



0.5 mm
Machined

Tapered 2°

Screw Pitch
1.2 mm

Tapered
6 mm



CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 8007	Ø 8.0	7.0 mm	Ø 5.0	Ø 6.0
IK 8008	Ø 8.0	8.5 mm	Ø 5.0	Ø 6.0
IK 8010	Ø 8.0	10 mm	Ø 5.0	Ø 6.0

ONLY FOR SOFT AND STD BONE - D3/D4

